

Stepwise Sweep VEP

Visual Acuity Assessment Module



The Diagnosis Stepwise Sweep VEP (ssVEP) module is an automated method to assess visual acuity objectively using visual evoked potential (VEP) techniques. The Diagnosis ssVEP module is based on methods developed and published by Bach, et al^{1,2,3} and runs on the Diagnosis E3 and Profile systems.

Features

Method

- Automatic noise-corrected response amplitudes
 - LaPlace channel and fourier transforms used
- Utilizes a pattern onset VEP methodology
- Utilizes a high-performance monitor
 - Same monitor used for pattern VEP and mfERG testing
- System and VEP techniques utilize standard ISCEV VEP methods

Process to use

- Protocol takes about 4 minutes per eye
- Easy to follow instructions, similar to ISCEV pattern VEP
- Visual acuity estimate and confidence limit

Applications

- Patients aged 5 or older, within a range of acuity from 20/20 to 20/400 (6/6 to 6/120)
- Unexplained visual loss
- Non-organic visual loss assessment
- Functional neurological symptom disorder
- Objective assessment when subjective is questionable



Sweep VEP protocol running on Diagnosis systems at 40% contrast



1. Bach, M, Maurer, JP, Wolf, ME; "Visual evoked potential-based acuity assessment in normal vision, artificially degraded vision, and in patients;" Br J Ophthalmol 92:396-403, 2008.
2. Hoffman, MB, Brands, J, Behrens-Baumann, W, Bach, M; "VEP-based acuity assessment in low vision;" Doc Ophthalmol 135:209-218, 2017.
3. Bach, M, Farmer, JD; "Evaluation of the 'Freiburg Acuity VEP' on Commercial Equipment;" Doc Ophthalmol (2020) 140:139-145

Visual Acuity Assessment Report

Stepwise Sweep VEP Report




Clear Eye Clinic, Lowell, MA

Patient: Hospital #: DOB: 1994 Age: 24 Years Gender: Test date: 8/3/2018

Right Eye (OD)

Snellen 20ft	97.5% Conf.	Calculation	Sfo (CPD)	Log(SFo)
20/22	> 20/45	Auto	15.71	1.20


Estimated visual acuity value and the 97.5% confidence level is automatically presented for each eye



Left Eye (OS)

Snellen 20ft	97.5% Conf.	Calculation	Sfo (CPD)	Log(SFo)
20/24	> 20/48	Auto	14.66	1.17

VEP amplitude data, automatically noise-corrected, plotted over the range of tested check sizes



Test comment: bare eye
 Test type: Sweep VEP Normal Vision 27in Rev B.2 (15777) Field of view: 30 wide x 18 high (degrees)
 Pupil diameter (DIOPTER) mm: 9.0 (3.0) Electrodes: Gold Cup Gold Cup
 Screen size: 150 (150) x 100 (100) cm: 1-0000 Electrodes: Gold Cup Gold Cup
 Printed: 5/12/2020 4:07:51 AM Page: 1 of 1

Stepwise Sweep VEP Module ordering information

- D370 Sweep VEP Add-on Module
- D370-1: for use on 27" monitors
- D370-2: for use on 32" monitors
- D370-3: for use on 23" monitors

Available as an add-on module for the following systems:

- Diagnosys E³ desktop systems
- Diagnosys Profile cart-based systems

Note: Diagnosys controller must be an ARM type, and Espion software must be v6.63 or later.

www.diagnosysllc.com

US: Diagnosys LLC; 55 Technology Drive, Suite 100, Lowell, MA 01851; 978-458-1600; sales@diagnosysllc.com
 EU: Diagnosys Vision Ltd; Office 117, DOC Building, Balheary Road, Swords, Dublin, K67 E5A0, Ireland; +44 (0) 1223 520699; mail@diagnosysvision.com
 UK: Diagnosys UK Ltd; 5 Trust Court, Chivers Way, Vision Park Histon, Cambridge, CB24 9PW, UK; +44 (0) 1223 520699; mail@diagnosysuk.co.uk

The information listed herein is accurate as of the date of printing, however may change at any time without notice. The contents may differ from the current status of approval of the product in your country. Please contact your local Diagnosys representative for more information. © 2021 by Diagnosys, LLC. All rights reserved.
 Doc: 15174 Rev: D ECN 1900 Date: 2 December 2022