

Clinical Systems

Objective Visual Function Testing

A Complete Suite of

Electroretinography (ERG)

Electroretinograms (ERGs) are recorded using corneal or skin electrodes to characterize dysfunction based on the size, shape, and speed of the retina's response

Full Field Electroretinogram (ffERG)

Evaluate rod and cone pathways

Inherited Retinal Diseases, Diabetic retinopathy


Pattern Electroretinogram (pERG)

Focal stimulation of macular cones and RGCs

Early Glaucoma detection, Macular Dystrophy

Normal full-field ERG


Dark and light adapted ERG steps illustrate rod and cone responses to flash stimuli of varying intensities.



Multifocal Electroretinogram (mfERG)


Topographical map of macular function

Macular Degeneration, Plaquenil Toxicity, Diabetic Retinopathy




Normal mfERG

Strong central peak with gradually sloping parafoveal response



Plaquenil Toxicity

Central peak preserved with significant loss of parafoveal function



Macular Degeneration

Severe central peak loss with residual parafoveal function


Electro-oculogram (EOG)

EOGs measure the standing potential between the cornea and RPE.

Electro-oculogram (EOG)

Retinal pigmented epithelium (RPE)

Bestophinopathy, AZOOR, AVMD



Normal EOG

A ratio of 2.0 between the RPE potential in the light as compared to in the dark

Visual Function Testing

Visually Evoked Potentials (VEP)

VEPs are performed with electrodes on the scalp above the visual cortex. Results can reveal delays or misrouting of the visual

Visually Evoked Potentials (VEP)


Detect disturbances in the visual pathway

Optic Neuropathy, Albinism, CVI

Sweep VEP (sVEP) Visual Acuity Assessment


Estimate visual acuity (VA) using pattern VEP methods

Functional Vision Loss, Unexplained Vision Loss, CVI



Normal pattern VEP

A robust positive peak at ~100 milliseconds (P100)



Normal flash VEP

A series of negative and positive waves useful in testing patients unable to perform a pattern VEP.

Psychophysical Tests

Psychophysical tests require patients to respond to visual stimuli with button boxes to assess light sensitivity.

Full-field Sensitivity Threshold (FST)


Quantify visual changes in low vision patients

Low Vision, IRD, Clinical Trials

Dark Adaptometry (DA)


Measure the time it takes to dark adapt

Macular Degeneration



High Quality FST

A light sensitivity threshold is determined when the patient reliably responds to a light stimulus



Dark Adaptometry Curve

Measures the time it takes to adapt from bright light to darkness

Diagnosys Systems

Diagnosys systems offer globally trusted, versatile diagnostic solutions with five plug-and-play stimulators and optional add-on modules.

Tabletop E3 Console



The **E3 Console** is a tabletop system that is the most widely used worldwide.

This system requires a table that is not included

Why Choose E3?

- Compact design for small spaces
- Portability option with travel case
- Cost-effective, globally trusted

Cart-based Profile



The **Profile** system cart can roll into any operating room, or be locked in place.

Why Choose Profile?

- Suitable for operating room testing
- Self-contained entirely on a cart
- Spatially flexible, movable on wheels

Diagnosys systems come in tabletop and cart-based variations. Both systems can interface with any stimulator or add-on module.



An example of a desktop system with 4 stimulators: the ColorDome, ColorFlash, Envoy, and LCD Monitor provide a complete testing suite



Two cart system examples, one LCD Monitor and one ColorDome on an operating room (OR) cart for testing supine or seated patients

Flash Stimulators



ColorDome

Mounted, binocular ganzfeld

- Infrared camera for patient monitoring in dark room
- Trillion-to-one luminance range that spans the entire range of human light perception



ColorFlash

Handheld, binocular, designed for pediatrics

- Testing can be done from a distance
- Attention-getting light and sound effects for children



ColorBurst

Handheld, monocular ganzfeld

- Well-suited for infants held in arms
- Ideal for Flash VEP

	ColorDome	ColorFlash	ColorBurst
Full-field ERG	●	●	●
PhNR	●	●	●
Flash VEP	●	●	●
Extended ISCEV	●	●	●
EOG	●	-	-
FST	● *	●	-
DA	● *	-	-
Pupillometry	● *	-	-

** Add-on module, not included with stimulator purchase*

Pattern Stimulators



LCD Monitor

Binocular LCD Monitor

- Flash artifact free with LI upgrade
- Fixation camera to monitor patient gaze



Envoy

Monocular, handheld OLED Monitor

- Artifact-free testing for easy PERG
- Can be used for bedside pattern testing

	LCD Monitor	Envoy
Multifocal ERG	●	-
Pattern ERG	● †	●
Pattern VEP	●	●
Sweep VEP (Visual Acuity)	● *,†	-

*† Requires LI-upgrade * Add-on module, not included with stimulator purchase*


Reliable, Precise, and Accurate

Our amplifiers detect the tiniest signals where other systems fail, enabling us to record ERGs on patients who would otherwise be presumed to have extinguished signals.



Ultra-Sensitive Amplifier


A fully differential, high-resolution 1 nanoVolt (nV), 32-bit amplifier with a 5 volt range and common mode rejection



Unmatched repeatability

Here are 44 overlaid PhNR responses recorded from the same person, across multiple Diagnosys systems in 8 different test sessions

Easily Interpretable Reports



Tests & Reference Data

- Systems include ISCEV standard protocols
- Normative reference data is available
- Easily edit and create custom protocols

Output Test Results

- Concise reports to simplify interpretation
- Customizable report layouts
- DICOM compatible outputs

Our ultra-sensitive amplifiers offer unmatched test repeatability

Our preferred corneal electrode is the **DTL Plus Electrode**; however, our systems work with any non-proprietary electrodes.



Fine conductive threads


- Well-tolerated
- Direct corneal contact
- Minimal risk of abrasion or infection

Customizable Espion Software




- ISCEV standard protocols are included; custom protocols can be created easily
- Software is pre-loaded with all systems; software updates are always free
- All systems are sold with all-in-one computers
- Encryption and audit trails ensure data security

Add-on Modules for expanded testing capabilities



Dark Adaptometry

An early biomarker of AMD, and an ideal replacement for the Goldmann Weekers Dark Adaptometer




DiagnosysFST

Quantify light perception in patients with very low vision once the ERG goes flat and as a clinical trial endpoint for interventional therapies



iMask

Restrict visual field size by creating a partial-field stimulus for any ColorDome test



Pupillometry

Measure the pupillary light reflex by recording changes in pupil diameter and area



Objective Visual Acuity (VA)

Quantify the estimated VA with established Sweep VEP methods for fixating and non-fixating patients and children



LI-LCD Monitor

Upgrade the LCD monitor to have isoluminant pattern reversal stimulation required for PERG and VA testing

The industry leader.

Decades of experience, thousands of installations.

- *Most experienced customer support team in the industry*
- *Hundreds of peer-reviewed papers published using our systems*
- *Extensive world-wide network of leading experts operating Diagnosys systems*

Pediatric Testing

*Diagnosys Exclusive
Abridged Pediatric Protocol*



- No formal dark adaptation
- No anesthesia
- Skin electrodes

Clinical Trials

Consultation Services

- Pre-study consultation
- Site setup and training
- Data Services



Ordering information

D310: Profile Cart-based System

D315: E3 Tabletop System

All systems are configured with selected stimulators and modules.

Please contact us if you would like more information on Diagnosys products.

US: Diagnosys LLC; 55 Technology Drive, Suite 100, Lowell, MA 01851; 978-458-1600; sales@diagnosysllc.com

EU: Diagnosys Vision Ltd; Office 117, DOC Building, Balheary Road, Swords, Dublin, K67 E5A0, Ireland; +44 (0) 1223 520699; mail@diagnosysvision.com

UK: Diagnosys UK Ltd; 5 Trust Court, Chivers Way, Vision Park Histon, Cambridge, CB24 9PW, UK; +44 (0) 1223 520699; mail@diagnosysuk.co.uk

The information listed herein is accurate as of the date of printing, however may change at any time without notice. The contents may differ from the current status of approval of the product in your country. Please contact your local Diagnosys representative for more information. © 2024 by Diagnosys, LLC. All rights reserved.